

Original Paper

The QardioArm App in the Assessment of Blood Pressure and Heart Rate: Reliability and Validity Study

Victoria Mazoterias Pardo^{1*}, RN; Marta E Losa Iglesias^{2*}, MSc, PhD; José López Chicharro^{1*}, MD, PhD; Ricardo Becerro de Bengoa Vallejo^{1*}, PhD

¹Universidad Complutense de Madrid, Madrid, Spain

²Universidad Rey Juan Carlos, Alcorcon, Spain

* all authors contributed equally

Corresponding Author:

Marta E Losa Iglesias, MSc, PhD

Universidad Rey Juan Carlos

Avda Atenas s/n

Alcorcon, 28922

Spain

Phone: 34 616962413

Email: marta.losa@urjc.es

Abstract

Background: Self-measurement of blood pressure is a priority strategy for managing blood pressure.

Objective: The aim of this study was to evaluate the reliability and validity of blood pressure and heart rate following the European Society of Hypertension's international validation protocol, as measured with the QardioArm, a fully automatic, noninvasive wireless blood pressure monitor and mobile app.

Methods: A total of 100 healthy volunteers older than 25 years from the general population of Ciudad Real, Spain, participated in a test-retest validation study with two measurement sessions separated by 5 to 7 days. In each measurement session, seven systolic blood pressure, diastolic blood pressure, and heart rate assessments were taken, alternating between the two devices. The test device was the QardioArm and the previously validated criterion device was the Omron M3. Sessions took place at a single study site with an evaluation room that was maintained at an appropriate temperature and kept free from noises and distractions.

Results: The QardioArm displayed very consistent readings both within and across sessions (intraclass correlation coefficients=0.80-0.95, standard errors of measurement=2.5-5.4). The QardioArm measurements corresponded closely to those from the criterion device ($r>.96$) and mean values for the two devices were nearly identical. The QardioArm easily passed all validation standards set by the European Society of Hypertension International Protocol.

Conclusions: The QardioArm mobile app has validity and it can be used free of major measurement error.

(*JMIR Mhealth Uhealth* 2017;5(12):e198) doi: [10.2196/mhealth.8458](https://doi.org/10.2196/mhealth.8458)

KEYWORDS

blood pressure; heart rate; reliability; validity; mobile apps

Introduction

Blood pressure (BP) plays a major role in the development of cardiovascular disease, which is the leading cause of premature death worldwide [1]. For patients who have or are at risk for high BP, regular monitoring of BP is useful [1,2].

Successful BP management depends to a large extent on the patient's willingness and capacity to make certain lifestyle

changes [3]. As demonstrated in studies by the UK National Health System, self-measurement of BP is a priority strategy for managing BP [4]. Self-measured BP monitoring has many advantages. It enables diagnosis of hypertension or hypotension, helps patients control BP, improves therapeutic compliance, and minimizes the "white coat" and "masked hypertension" syndromes and observation biases caused by the health care professional being aware of the patient's characteristics during the measurement [5-8].

Figure 1. QardioArm is a device with a cuff used with an app for visualization of blood pressure known as Qardio.



Therefore, increased use of self-measured BP in the care of patients with arterial hypotension or hypertension could be beneficial. Recent clinical practice guidelines identified several indications for self-measured BP, which may provide for more standardized measurement of BP [5-17]. Today's powerful mobile technologies can greatly facilitate self-measured BP by giving individuals more control over their health and well-being than conventional clinic-bound medical devices and allowing clinicians to evaluate the patient's status and provide timely feedback remotely [18-28].

In a recent market study of mobile phone apps for managing hypertension, most were focused on health management. Only 14% of these Android and iPhone apps transformed a mobile phone into a BP measuring device [22]. QardioArm (Atten Electronic Co Ltd, Dongguan, China) is a device with a cuff used with an app for visualization known as Qardio [29] (Figure 1), and it may be a good alternative to enable self-measurement of BP by patients. None of these apps involved a BP cuff nor had any been validated against a gold standard. Only 3% of these apps were developed by health organizations, such as universities or professional organizations [18,20,22,23,25-28]. Nonetheless, consumers downloaded these apps for measuring BP and heart rate (HR) and evaluated them favorably, even though they had not been validated [22,24,26,30].

Self-measured BP devices are only useful and beneficial to the extent that they are user-friendly and accurate [6-9,11,13,14]. Self-measured BP devices should be validated by independent experts according to accepted protocols designed specifically for this purpose, such as those established by the British Hypertension Society [31], the Association for the Advancement of Medical Instrumentation [32], and the European Society of Hypertension (ESH) [33,34]. The ESH called for self-measured BP devices to be validated with the 2010 version of the ESH

International Protocol [35]. Accordingly, we evaluated the accuracy of self-measured BP and HR measurements obtained with one such mobile app, the QardioArm, with the ESH International Protocol.

Methods

We examined the concordance between two European Community-certified devices for measuring BP and HR in a study involving repeated measurements conducted in the city of Ciudad Real, Spain, between May and August 2016. The Clinical Research Ethics Committee of Hospital Clínico San Carlos in Madrid, Spain, approved this study (number 16/179-E). This study complies with the ethical principles of the Declaration of Helsinki [36], including amendments from 2000 to 2013.

Participants

We recruited 125 participants from the general population of Ciudad Real using snowball sampling, from May to August 2016. A total of 100 participants completed the study; 25 participants did not return for the second measurement. According to the ESH International Protocol [35], a minimum of one-third of the total number of participants must be men and a minimum of one-third must be women. The participants satisfied the eligibility criteria in the latest revision of the ESH International Protocol [33-35], which are as follows:

Inclusion criteria:

1. Demographic: men and women aged at least 25 years. Of the total number of participants, at least one-third must be men and one-third must be women.
2. Temporal and geographic: healthy adult volunteers able to attend assessments in Ciudad Real.
3. Clinical: ideal BP range of 90 to 180 mm Hg for systolic blood pressure (SBP) and 40 to 130 mm Hg for diastolic

blood pressure (DBP), and arm circumference of 220 to 320 mm.

Exclusion criteria:

1. Clinical: sustained arrhythmia, circulatory problems in which the use of a cuff is contraindicated, and/or pregnancy.
2. Cognitive impairment that leads to inability to follow instructions.

All participants gave their written informed consent to participate in this study.

Before recruitment, we calculated with the GRANMO sample size calculation program (Institut Hospital del Mar d'Investigacions Mèdiques, Barcelona, Spain) that a sample size of 125 would give 80% statistical power ($\alpha=.05$) and two-tailed test to detect a difference equal to or greater than 0.9 mm Hg in BP. This calculation was based on the assumptions of a standard deviation of 3.2 mm Hg [37] and a loss to follow-up rate of 20%.

Procedure

Participants visited the study site for two measurement sessions separated by 5 to 7 days. The first author (VMP) collected all the data. Both BP and HR were measured in the left arm. The central portion of the cuff was placed at heart level, with the bottom edge 1 to 2 cm above the elbow.

During the first session, each participant reported his/her sex and date of birth. Measurements were taken of the participant's weight (with an automatic digital scale), height (with a wall-mounted measuring rod), and arm circumference (with a tape measure at the midpoint between the acromion and olecranon). The participant then relaxed for 10 minutes before baseline BP and HR were measured with the criterion BP

monitor Omron M3 (HEM-7200-E, Omron Healthcare Co Kyoto, Japan) (Figure 2). Next, BP and HR were measured with the test device QardioArm only to confirm that the test device was working correctly on the participant (these were not included in analysis).

The rest of the first session involved seven BP and HR measurements, alternating between the criterion and test devices, with 2 to 3 minute intervals between measurements. The first, third, fifth, and seventh measurements were made with the criterion device; the second, fourth, and sixth measurements were made with the test device. The second measurement session involved the same series of seven measurements, with identical procedures carried out by the same researcher, which occurred on the same day of the week and at the same time after lunch, in the same room with identical conditions, as far as it was possible.

During measurements, participants were calm and quiet while sitting with their feet parallel and flat on the floor, their legs uncrossed, and their left hands resting palm side up on a flat surface. The room was maintained at an appropriate temperature and kept free from noises and distractions during the sessions [33,34].

Of the 125 participants who completed the first session, 100 also completed the second session. Our analysis focused on data from just these 100 participants.

Study Devices

The Omron M3 Intellisense was the criterion device used as a benchmark in our study. It has CE 0197 certification and has been validated by the ESH between the Omron M3 and a different criterion device with mean 1.7 (SD 3.2) mm Hg for systolic and mean -0.9 (SD 2.6) mm Hg for diastolic [36].

Figure 2. The Omron M3 blood pressure monitor.



The Omron M3 is a fully automatic BP monitor operating on the oscillometric principle. It has a range of 0 to 299 mm Hg for BP and 40 to 180 beats per minute (bpm) for HR. The cuff is inflated using an electric pump and deflated by means of a pressure release valve.

After each measurement, the SBP, DBP, and HR are shown on the Omron M3's LCD screen. The device can also display a symbol on the screen indicating an irregular heartbeat, which is detected during measurement of SBP and DBP.

The QardioArm app for mobile phones and tablets was the test device in our study. It has CE 0734 certification and is a fully automatic, noninvasive wireless BP monitor for measuring SBP, DBP, and HR in adults. The QardioArm also operates on the oscillometric principle and has an inflatable cuff, which is placed around the upper arm. The cuff is inflated automatically and deflated using a controlled pressure release valve. The cuff is suitable for an arm circumference of between 22 and 37 cm. The measuring range is 40 to 250 mm Hg for BP and 40 to 200 bpm for HR.

The QardioArm has an automatic screen with graphics and visuals to facilitate interpretation of the data. The monitor connects to the free Qardio app on any device with Bluetooth 4.0 that runs iOS 7.0 or later or Android 4.4 or later. It can be used with iPhones, iPods, iPads, Apple Watches, and Android mobile phones and tablets. The data are stored on the mobile phone or tablet and therefore the pattern of values over time can be viewed. The app can be configured to issue reminders and warnings, and the measurements and progress can be shared in real time with other users.

The QardioArm weighs approximately 310 g without batteries and measures 6.8 cm in width, 3.8 cm in height, and 14 cm in length, with the cuff closed. It requires four AAA batteries.

Statistical Analysis

All data were analyzed using SPSS version 19 (IBM Corp, Armonk, NY, USA). Statistical significance was set at $P < .05$.

We computed univariate summary statistics on participants' characteristics. Body mass index (BMI) was defined as weight in kilograms divided by height (in meters) squared. To assess whether the BP and HR variables were normally distributed, we conducted one-sample Kolmogorov-Smirnov tests, applying the Lilliefors significance correction. Significant probability values in this test indicate nonnormality.

To evaluate the reliability of the QardioArm, we computed the intraclass correlation coefficient (ICC) for comparing measurements taken within the same session (intrasession) and in different sessions (intersession) [38]. The participants were the groups in each ICC analysis. Each intrasession evaluation involved the three measurements obtained using the test device for a particular primary variable (SBP, DBP, and HR). We used the session means of a variable in calculating the intersession ICC.

We interpreted ICC values with the guide proposed by Landis and Koch [39]: ≤ 0.20 =slight agreement, $0.21-0.40$ =fair agreement, $0.41-0.60$ =moderate agreement, $0.61-0.80$ =substantial agreement, and ≥ 0.81 =almost perfect

agreement. We also followed the recommendations of Portney and Watkins [40], who suggested that an ICC of 0.90 or higher for clinical measurements indicates they are reliable.

We also calculated the standard error of measurement (SEM) for each intrasession and intersession primary variable. SEM values indicate the degree of error, with low values indicating low error. We used the SEM formula of the standard deviation multiplied by the square root of $(1-ICC)$. To identify systematic error over time, we also compared first and second session means on each primary variable with paired-sample Student t tests.

To determine the validity of the QardioArm, we compared the QardioArm measurements with those made by the Omron M3, using both the ESH's International Protocol for the validation of BP measuring devices and the Student t test. Following the ESH protocol [33,34], we computed the absolute values of the differences between successive pairs of the seven measurements in a session (second–first, second–third, fourth–third, fourth–fifth, sixth–fifth, and sixth–seventh). This gives three paired differences for each variable (SBP, DBP, and HR), participant, and session. We classified whether the paired differences were ≤ 5 , ≤ 10 , ≤ 15 , or > 15 mm Hg for BP and ≤ 3 , ≤ 5 , ≤ 8 , or > 8 bpm in the case of HR. The four levels of difference correspond to very accurate, slightly inaccurate, moderately inaccurate, and very inaccurate, respectively [33,34]. If the device passed during both sessions, it could be recommended for clinical use. If it did not pass during the second session, it could not be recommended for clinical use.

We compared the observed number of differences (out of 300: three paired differences for a variable \times 100 participants in a session) falling into these categories with the standards specified in the ESH International Protocol. There were two phases of this comparison. If a test device did not pass during the first session, the validation process was terminated. If it passed during the first session, the second session was started.

In addition to following the ESH International Protocol, we computed Student t test for independent samples to check whether there were differences between devices in BP and HR measurements. For these tests, we used the mean values of the three study variables obtained from both sessions for the different devices. We also compared the mean intersession differences for the two devices. Furthermore, we calculated the Pearson correlations between the devices for the mean participant values on the three primary variables [41].

Finally, we produced Bland-Altman plots [42] to display the agreement between the two devices. These plots show the difference between each pair of measurements on the y-axis against the mean of each pair of measurements on the x-axis.

Results

Table 1 shows the sociodemographic characteristics of the study participants ($N=100$). There were significant mean differences between men and women in weight ($P < .001$), height ($P < .001$), and BMI ($P = .009$), but no differences in age ($P = .10$) and arm circumference ($P = .08$). The Kolmogorov-Smirnov tests showed

that SBP ($P=.14$), DBP ($P=.20$), HR ($P=.20$), as well as weight ($P=.17$) and BMI ($P=.20$) were normally distributed

Table 2 shows the intrasession and intersession reliability results for the QardioArm. In every case, reliability was high for SBP, DBP, and HR, and there was no evidence of systematic error from one session to the next.

Table 3 shows the validity results for the QardioArm, with the Omron M3 as the criterion. The QardioArm passed both sessions of the ESH validation process for SBP, DBP, and HR.

Table 4 shows the comparisons between the devices in terms of mean values on the primary variables as well as the Pearson

correlations between them at the participant level. The two devices produced very similar mean BP and HR estimates. The differences in mean BP values for the QardioArm and Omron M3 were less than 2 mm Hg and the differences in mean HR values between devices were less than 0.1 bpm. All intrasession and intersession Pearson correlations were greater than .96.

Figures 3, 4, and 5 display the Brand-Altman plots for the BP and HR variables in the first session. For each variable and almost every participant, the difference between device means fell within the 95% confidence interval of all measurements. The plots for session 2 (not shown) gave very similar results.

Table 1. Sociodemographic and clinical characteristics of the sample.

Variables	Total group (N=100)		Male (n=37)		Female (n=63)	
	Mean (SD)	Range	Mean (SD)	Range	Mean (SD)	Range
Age (years)	52.2 (18.8)	25.0-89.0	52.3 (20.9)	25.0-89.0	52.2 (17.6)	25.0-83.0
Weight (kg)	76.2 (14.6)	49.5-125.1	84.6 (13.1)	65.3-125.1	71.3 (13.2)	49.5-102.4
Height (cm)	167.3 (8.1)	151.0-187.0	172.1 (7.1)	160.0-184.0	164.5 (7.4)	151.0-187.0
BMI ^a (kg/m ²)	27.1 (4.2)	18.9-38.2	28.5 (3.8)	22.2-38.2	26.3 (4.2)	18.9-37.7
Arm circumference (mm)	295.5 (35.2)	234.0-320.0	302.7 (33.9)	245.0-320.0	291.4 (35.6)	234.0-320.0
Baseline SBP^b						
Omron	130.9 (17.3)	88.0-180.0	136.0 (17.7)	101.0-180.0	127.9 (16.4)	88.0-167.0
QardioArm	128.3 (17.2)	87.0-178.0	131.9 (17.4)	100.0-178.0	126.1 (16.8)	87.0-168.0
Baseline DBP^c						
Omron	74.2 (10.0)	49.0-98.0	73.5 (11.7)	57.0-98.0	74.6 (8.9)	49.0-97.0
QardioArm	75.8 (10.8)	49.0-103.0	75.4 (12.3)	54.0-103.0	75.9 (9.9)	49.0-103.0

^aBMI: body mass index.

^bSBP: systolic blood pressure.

^cDBP: diastolic blood pressure.

Table 2. Intrasession and intersession reliability for the QardioArm in sessions 1 and 2.

Variable	Mean (SD)	Range	ICC (95% CI)	P value	SEM
Session 1					
Systolic pressure (mm Hg)	126.5 (16.1)	87.0-179.0	0.89 (0.84-0.92)		5.35
Diastolic pressure (mm Hg)	74.8 (10.3)	46.7-100.3	0.91 (0.87-0.93)		3.15
Heart rate (bpm)	71.4 (10.4)	53.3-100.7	0.92 (0.87-0.94)		3.00
Session 2					
Systolic pressure (mm Hg)	124.9 (15.8)	89.7-171.7	0.91 (0.86-0.95)		4.67
Diastolic pressure (mm Hg)	73.9 (9.6)	51.7-95.7	0.87 (0.82-0.90)		3.52
Heart rate (bpm)	70.1 (11.0)	44.3-96.0	0.95 (0.93-0.96)		2.53
Intersession					
Systolic pressure (mm Hg)	125.7 (15.3)	87.0-179.0	0.91 (0.87-0.94)	.09	4.73
Diastolic pressure (mm Hg)	74.3 (9.4)	46.7-100.3	(0.83-0.92)	.17	3.36
Heart rate (bpm)	70.7 (9.8)	44.3-100.7	0.80 (0.70-0.86)	.14	4.82

Table 3. shows the validity results for the QardioArm, with the Omron M3 as the criterion. The QardioArm passed both sessions of the ESH validation process for SBP, DBP, and HR.

	Category					Result	Mean difference \pm SD ^a
	5 mm Hg/ 3 bpm ^b	10 mm Hg/ 5 bpm ^b	15 mm Hg/ 8 bpm ^b	2/3 within 5 mm Hg/3 bpm ^c	0/3 within 5 mm Hg/3 bpm ^c		
	Session 1						
Required							
Two of ^d	219	261	288	N/A ^e	N/A	N/A	N/A
All of ^f	195	243	279	N/A	N/A	N/A	N/A
Observed							
SBP1 ^{g,h}	242	286	298	N/A	N/A	Pass	3.43 \pm 3.13
DBP1 ⁱ	277	295	299	N/A	N/A	Pass	2.43 \pm 2.35
HR1 ^j	265	287	297	N/A	N/A	Pass	1.60 \pm 1.92
SBP2 ^k	256	289	300	N/A	N/A	Pass	3.29 \pm 2.74
DBP2	266	295	298	N/A	N/A	Pass	2.79 \pm 2.86
HR2	272	287	295	N/A	N/A	Pass	1.61 \pm 1.99
Session 2							
Required	N/A	N/A	N/A	≥ 72 ^l	≤ 9 ^m	N/A	N/A
Observed							
SBP1	N/A	N/A	N/A	87	0	Pass	N/A
DBP1	N/A	N/A	N/A	95	0	Pass	N/A
HR1	N/A	N/A	N/A	90	0	Pass	N/A
SBP2	N/A	N/A	N/A	94	0	Pass	N/A
DBP2	N/A	N/A	N/A	95	2	Pass	N/A
HR2	N/A	N/A	N/A	93	0	Pass	N/A

^aSD: standard deviation.

^bFrequencies of paired measurements (of 300) within particular ranges.

^cFrequencies of participants (of 100) with measurements meeting certain criteria.

^d“Two of” means that two of the three conditions given in the same row must be met.

^eN/A: not applicable.

^f“All of” means that all three requirements given in the row must be met.

^gSBP: systolic blood pressure.

^hSession 1.

ⁱDBP: diastolic blood pressure.

^jHR: heart rate.

^kSession 2.

^lAt least two of the three pairs of differences for at least 72 participants must fall in the 5 mmHg category.

^mThe three pairs of differences of no more than 9 participants may fall outside the 5 mmHg category.

Table 4. Comparison of Qardioarm and Omron devices in mean estimates and correlations of individual means between devices.

Variable	Omron, mean (SD) ^a	QardioArm, mean (SD)	t-test <i>P</i> value	<i>r</i>	<i>P</i> value
Intrasession 1					
Systolic pressure (mm Hg)	127.34 (16.02)	126.47 (16.06)	.70	.983	<.001
Diastolic pressure (mm Hg)	73.66 (9.82)	74.78 (10.29)	.43	.984	<.001
Heart rate (bpm)	71.31 (10.29)	71.36 (10.42)	.97	.987	<.001
Intrasession 2					
Systolic pressure (mm Hg)	125.65 (16.02)	124.96 (15.84)	.76	.984	<.001
Diastolic pressure (mm Hg)	72.23 (9.08)	73.90 (9.61)	.21	.966	<.001
Heart rate (bpm)	70.12 (10.65)	70.05 (10.99)	.96	.987	<.001
Intersession					
Systolic pressure (mm Hg)	126.50 (15.47)	125.71 (15.30)	.72	.991	<.001
Diastolic pressure (mm Hg)	72.95 (8.90)	74.34 (9.43)	.29	.982	<.001
Heart rate (bpm)	70.72 (9.68)	70.70 (9.76)	.99	.993	<.001

^aSD: standard deviation.


Figure 3. Bland-Altman plot comparing QardioArm and Omron devices for systolic blood pressure (SBP), for individual participants in session 1.

Figure 4. Bland-Altman plot comparing QardioArm and Omron devices for diastolic blood pressure (DBP) for individual participants in session 1.





Figure 5. Bland-Altman plot comparing QardioArm and Omron devices for heart rate (HR) for individual participants in session 1.



Discussion

The QardioArm is recommended by dabl Educational Trust, which is an organization that provides an interpreted one-page BP report for the medical team and patient [43], but there are not any validation studies on the QardioArm published in a peer-reviewed journal. Recently, a study compared QardioArm with three other mobile phone-compatible BP measuring devices using a handheld aneroid sphygmomanometer as the reference device [44], but the study did not follow the ESH International Protocol.

We evaluated the reliability and validity of the QardioArm device for self-measured BP with 100 participants in two

measurement sessions. The QardioArm displayed very consistent readings both within and across sessions. The QardioArm measurements corresponded closely to those from the previously validated criterion device, the Omron M3. The QardioArm easily passed all validation standards set by the ESH International Protocol. Therefore, the QardioArm can be recommended for clinical use in individuals with similar characteristics to those who participated in this study, such as adults aged 48 to 56 years with a BMI between 25 and 29, not pregnant, and without cardiac arrhythmia, vascular problems, or arteriovenous fistulas in the arm.

There are several limitations of the ESH International Protocol and, thus, our study. The ESH International Protocol uses the

same validation requirements (eg, ≤ 5 mm Hg difference indicating accuracy) for both SBP and DBP, even though the magnitude of diastolic values is usually half that of systolic measurements.

Also, HR as a parameter is not considered in any version of the ESH International Protocol despite the fact that automatic sphygmomanometers produce measurements of HR as well as BP. We have not found any studies of HR validation attempts in the literature that follow the ESH International Protocol [37,45-49]. Therefore, ours may be the first study to validate this parameter with the criteria of the ESH protocol for BP. We set the requirements for HR validation to be roughly proportional to the magnitude of the measurements, and they may even be more demanding than the BP requirements of the ESH International Protocol.

Moreover, the recommendations made in the ESH International Protocol regarding the populations to which validation results can be applied may not be followed faithfully in clinical practice. The ESH International Protocol imposes certain gender requirements and limits validation studies to individuals older than 25 years who have BPs and arm circumferences within specific ranges. Because these subgroups represent only a portion of the large heterogeneous population with BP abnormalities, extrapolation of ESH validation results to other specific populations could be considered unsafe. Additional

validation studies are needed if the product is to be used in other subpopulations, such as pregnant women, obese individuals, children, or individuals with particular conditions, such as arrhythmia.

The ESH International Protocol also does not indicate the number of validation studies needed to confirm the accuracy of the device. According to experts, at least two validation studies should be performed in different centers and in different populations [37]. The protocol of the Association for the Advancement of Medical Instrumentation recommends more than one study, but does not specify the number of studies or devices [32]. Because the QardioArm device has not been validated previously, we are unable to compare our findings with those of other authors. When compared with the Omron M3, the QardioArm tends to give very slightly and nonsignificantly lower SBP and very slightly and nonsignificantly higher DBP readings. We recommend that the QardioArm be further validated with different study designs and study sites and with different types of populations.

There is a high concordance between the measurements made with QardioArm, both intrasession and intersession. The QardioArm mobile app has validity and there is a direct linear correlation between QardioArm measurements and the previously validated Omron M3 measurements, so it can be used without the risk of major measurement error.

Conflicts of Interest

None declared.

References

1. Wang H, Dwyer-Lindgren L, Lofgren KT, Rajaratnam JK, Marcus JR, Levin-Rector A, et al. Age-specific and sex-specific mortality in 187 countries, 1970–2010: a systematic analysis for the Global Burden of Disease Study 2010. *Lancet* 2012 Dec;380(9859):2071-2094. [doi: [10.1016/S0140-6736\(12\)61719-X](https://doi.org/10.1016/S0140-6736(12)61719-X)]
2. Chow CK, Teo KK, Rangarajan S, Islam S, Gupta R, Avezum A, PURE (Prospective Urban Rural Epidemiology) Study investigators. Prevalence, awareness, treatment, and control of hypertension in rural and urban communities in high-, middle-, and low-income countries. *JAMA* 2013 Sep 04;310(9):959-968. [doi: [10.1001/jama.2013.184182](https://doi.org/10.1001/jama.2013.184182)] [Medline: [24002282](https://pubmed.ncbi.nlm.nih.gov/24002282/)]
3. Khatib R, Schwalm J, Yusuf S, Haynes RB, McKee M, Khan M, et al. Patient and healthcare provider barriers to hypertension awareness, treatment and follow up: a systematic review and meta-analysis of qualitative and quantitative studies. *PLoS One* 2014;9(1):e84238 [FREE Full text] [doi: [10.1371/journal.pone.0084238](https://doi.org/10.1371/journal.pone.0084238)] [Medline: [24454721](https://pubmed.ncbi.nlm.nih.gov/24454721/)]
4. Department of Health. The NHS Outcomes Framework 2012/13. London: Department of Health; Dec 09, 2011:13.
5. McManus RJ, Mant J, Haque MS, Bray EP, Bryan S, Greenfield SM, et al. Effect of self-monitoring and medication self-titration on systolic blood pressure in hypertensive patients at high risk of cardiovascular disease: the TASMIN-SR randomized clinical trial. *JAMA* 2014 Aug 27;312(8):799-808. [doi: [10.1001/jama.2014.10057](https://doi.org/10.1001/jama.2014.10057)] [Medline: [25157723](https://pubmed.ncbi.nlm.nih.gov/25157723/)]
6. Uhlig K, Patel K, Ip S, Kitsios GD, Balk EM. Self-measured blood pressure monitoring in the management of hypertension: a systematic review and meta-analysis. *Ann Intern Med* 2013 Aug 06;159(3):185-194. [doi: [10.7326/0003-4819-159-3-201308060-00008](https://doi.org/10.7326/0003-4819-159-3-201308060-00008)] [Medline: [23922064](https://pubmed.ncbi.nlm.nih.gov/23922064/)]
7. Sharman JE, Howes FS, Head GA, McGrath BP, Stowasser M, Schlaich M, et al. Home blood pressure monitoring: Australian Expert Consensus Statement. *J Hypertens* 2015 Sep;33(9):1721-1728 [FREE Full text] [doi: [10.1097/HJH.0000000000000673](https://doi.org/10.1097/HJH.0000000000000673)] [Medline: [26136205](https://pubmed.ncbi.nlm.nih.gov/26136205/)]
8. Howes F, Hansen E, Nelson M. Management of hypertension in general practice--a qualitative needs assessment of Australian GPs. *Aust Fam Physician* 2012 May;41(5):317-323 [FREE Full text] [Medline: [22558624](https://pubmed.ncbi.nlm.nih.gov/22558624/)]
9. Parati G, Stergiou GS, Asmar R, Bilo G, de LP, Imai Y, ESH Working Group on Blood Pressure Monitoring. European Society of Hypertension practice guidelines for home blood pressure monitoring. *J Hum Hypertens* 2010 Dec;24(12):779-785. [doi: [10.1038/jhh.2010.54](https://doi.org/10.1038/jhh.2010.54)] [Medline: [20520631](https://pubmed.ncbi.nlm.nih.gov/20520631/)]

10. Niiranen TJ, Johansson JK, Reunanen A, Jula AM. Optimal schedule for home blood pressure measurement based on prognostic data: the Finn-Home Study. *Hypertension* 2011 Jun;57(6):1081-1086 [FREE Full text] [doi: [10.1161/HYPERTENSIONAHA.110.162123](https://doi.org/10.1161/HYPERTENSIONAHA.110.162123)] [Medline: [21482956](https://pubmed.ncbi.nlm.nih.gov/21482956/)]
11. Mengden T, Dechend R, Lüders S, Eckert S, Tholl U, Sanner B. [Practice guidelines for home blood pressure measurement]. *Dtsch Med Wochenschr* 2011 Nov;136(46):2374-2377. [doi: [10.1055/s-0031-1292056](https://doi.org/10.1055/s-0031-1292056)] [Medline: [22068451](https://pubmed.ncbi.nlm.nih.gov/22068451/)]
12. Daskalopoulou S, Khan N, Quinn R, Ruzicka M, McKay D, Hackam D, Canadian Hypertension Education Program. The 2012 Canadian hypertension education program recommendations for the management of hypertension: blood pressure measurement, diagnosis, assessment of risk, and therapy. *Can J Cardiol* 2012 May;28(3):270-287. [doi: [10.1016/j.cjca.2012.02.018](https://doi.org/10.1016/j.cjca.2012.02.018)] [Medline: [22595447](https://pubmed.ncbi.nlm.nih.gov/22595447/)]
13. Imai Y, Kario K, Shimada K, Kawano Y, Hasebe N, Matsuura H, Japanese Society of Hypertension Committee for Guidelines for Self-monitoring of Blood Pressure at Home. The Japanese Society of Hypertension Guidelines for Self-monitoring of Blood Pressure at Home (Second Edition). *Hypertens Res* 2012 Aug;35(8):777-795. [doi: [10.1038/hr.2012.56](https://doi.org/10.1038/hr.2012.56)] [Medline: [22863910](https://pubmed.ncbi.nlm.nih.gov/22863910/)]
14. Mancia G, Fagard R, Narkiewicz K, Redón J, Zanchetti A, Böhm M. ESH/ESC Guidelines for the management of arterial hypertension: The Task Force for the management of arterial hypertension of the European Society of Hypertension (ESH) and of the European Society of Cardiology (ESC). *J Hypertens* 2013;31(10):1281-1237. [doi: [10.1097/01.hjh.0000431740.32696.cc](https://doi.org/10.1097/01.hjh.0000431740.32696.cc)]
15. Boivin J, Boutte E, Fay R, Rossignol P, Zannad F. Home blood pressure monitoring: a few minutes of rest before measurement may not be appropriate. *Am J Hypertens* 2014 Jul;27(7):932-938. [doi: [10.1093/ajh/hpu001](https://doi.org/10.1093/ajh/hpu001)] [Medline: [24561656](https://pubmed.ncbi.nlm.nih.gov/24561656/)]
16. Niiranen TJ, Asayama K, Thijs L, Johansson JK, Hara A, Hozawa A, IDHOCO Investigators. Optimal number of days for home blood pressure measurement. *Am J Hypertens* 2015 May;28(5):595-603. [doi: [10.1093/ajh/hpu216](https://doi.org/10.1093/ajh/hpu216)] [Medline: [25399016](https://pubmed.ncbi.nlm.nih.gov/25399016/)]
17. Sharman JE, Howes F, Head GA, McGrath BP, Stowasser M, Schlaich M, et al. How to measure home blood pressure: recommendations for healthcare professionals and patients. *Aust Fam Physician* 2016;45(1):31-34 [FREE Full text] [Medline: [27051984](https://pubmed.ncbi.nlm.nih.gov/27051984/)]
18. de Jongh T, Gurol-Urganci I, Vodopivec-Jamsek V, Car J, Atun R. Mobile phone messaging for facilitating self-management of long-term illnesses. *Cochrane Database Syst Rev* 2012 Dec 12;12:CD007459. [doi: [10.1002/14651858.CD007459.pub2](https://doi.org/10.1002/14651858.CD007459.pub2)] [Medline: [23235644](https://pubmed.ncbi.nlm.nih.gov/23235644/)]
19. Mann D, Kudesia V, Reddy S, Weng M, Imler D, Quintiliani L. Development of DASH Mobile: a mHealth lifestyle change intervention for the management of hypertension. *Stud Health Technol Inform* 2013;192:973. [Medline: [23920747](https://pubmed.ncbi.nlm.nih.gov/23920747/)]
20. Martínez-Pérez B, de la Torre-Díez I, López-Coronado M, Herreros-González J. Mobile apps in cardiology: review. *JMIR Mhealth Uhealth* 2013 Jul 24;1(2):e15 [FREE Full text] [doi: [10.2196/mhealth.2737](https://doi.org/10.2196/mhealth.2737)] [Medline: [25098320](https://pubmed.ncbi.nlm.nih.gov/25098320/)]
21. Bradbury K, Watts S, Arden-Close E, Yardley L, Lewith G. Developing digital interventions: a methodological guide. *Evid Based Complement Alternat Med* 2014;2014:561320 [FREE Full text] [doi: [10.1155/2014/561320](https://doi.org/10.1155/2014/561320)] [Medline: [24648848](https://pubmed.ncbi.nlm.nih.gov/24648848/)]
22. Kumar N, Khunger M, Gupta A, Garg N. A content analysis of smartphone-based applications for hypertension management. *J Am Soc Hypertens* 2015 Feb;9(2):130-136. [doi: [10.1016/j.jash.2014.12.001](https://doi.org/10.1016/j.jash.2014.12.001)] [Medline: [25660364](https://pubmed.ncbi.nlm.nih.gov/25660364/)]
23. McLean G, Murray E, Band R, Saunderson K, Hanlon P, Little P, et al. Digital interventions to promote self-management in adults with hypertension: protocol for systematic review and meta-analysis. *JMIR Res Protoc* 2015 Nov 20;4(4):e133 [FREE Full text] [doi: [10.2196/resprot.4648](https://doi.org/10.2196/resprot.4648)] [Medline: [26589728](https://pubmed.ncbi.nlm.nih.gov/26589728/)]
24. Hallberg I, Ranerup A, Kjellgren K. Supporting the self-management of hypertension: patients' experiences of using a mobile phone-based system. *J Hum Hypertens* 2016 Feb;30(2):141-146 [FREE Full text] [doi: [10.1038/jhh.2015.37](https://doi.org/10.1038/jhh.2015.37)] [Medline: [25903164](https://pubmed.ncbi.nlm.nih.gov/25903164/)]
25. Omboni S, Caserini M, Coronetti C. Telemedicine and mHealth in hypertension management: technologies, applications and clinical evidence. *High Blood Press Cardiovasc Prev* 2016 Sep;23(3):187-196. [doi: [10.1007/s40292-016-0143-6](https://doi.org/10.1007/s40292-016-0143-6)] [Medline: [27072129](https://pubmed.ncbi.nlm.nih.gov/27072129/)]
26. Plante TB, Urrea B, MacFarlane ZT, Blumenthal RS, Miller ER, Appel LJ, et al. Validation of the instant blood pressure smartphone app. *JAMA Intern Med* 2016 May 01;176(5):700-702 [FREE Full text] [doi: [10.1001/jamainternmed.2016.0157](https://doi.org/10.1001/jamainternmed.2016.0157)] [Medline: [26938174](https://pubmed.ncbi.nlm.nih.gov/26938174/)]
27. Kang H, Park H. A mobile app for hypertension management based on clinical practice guidelines: development and deployment. *JMIR Mhealth Uhealth* 2016 Feb 02;4(1):e12 [FREE Full text] [doi: [10.2196/mhealth.4966](https://doi.org/10.2196/mhealth.4966)] [Medline: [26839283](https://pubmed.ncbi.nlm.nih.gov/26839283/)]
28. İlhan I, Yıldız I, Kayrak M. Development of a wireless blood pressure measuring device with smart mobile device. *Comput Methods Programs Biomed* 2016 Mar;125:94-102. [doi: [10.1016/j.cmpb.2015.11.003](https://doi.org/10.1016/j.cmpb.2015.11.003)] [Medline: [26626086](https://pubmed.ncbi.nlm.nih.gov/26626086/)]
29. Qardio. 2013. URL: <https://www.getqardio.com/es/> [accessed 2017-11-30] [WebCite Cache ID 6vMQzFn83]
30. dabl Educational Trust. 2015. Devices for blood pressure measurement URL: <http://www.dableducational.org/> [accessed 2017-12-04] [WebCite Cache ID 6vSem9Wkl]
31. O'Brien E, Petrie J, Littler W, de SM, Padfield P, Altman D, et al. An outline of the revised British Hypertension Society protocol for the evaluation of blood pressure measuring devices. *J Hypertens* 1993 Jun;11(6):677-679. [Medline: [8397248](https://pubmed.ncbi.nlm.nih.gov/8397248/)]

32. American National Standard. ANSI/AAMI/ ISO 81060-2. Non-Invasive Sphygmomanometers. Part 2: Clinical Validation of Automated Measurement Type. Arlington, VA: Association for the Advancement of Medical Instrumentation URL: http://my.aami.org/aamiresources/previewfiles/8106002_1306_preview.pdf [accessed 2017-12-04] [WebCite Cache ID 6vSfQ23n2]
33. O'Brien E, Pickering T, Asmar R, Myers M, Parati G, Staessen J, Working Group on Blood Pressure Monitoring of the European Society of Hypertension. Working Group on Blood Pressure Monitoring of the European Society of Hypertension International Protocol for validation of blood pressure measuring devices in adults. *Blood Press Monit* 2002 Feb;7(1):3-17. [Medline: 12040236]
34. O'Brien E, Atkins N, Stergiou G, Karpettas N, Parati G, Asmar R, Working Group on Blood Pressure Monitoring of the European Society of Hypertension. European Society of Hypertension International Protocol revision 2010 for the validation of blood pressure measuring devices in adults. *Blood Press Monit* 2010 Feb;15(1):23-38. [doi: 10.1097/MBP.0b013e3283360e98] [Medline: 20110786]
35. O'Brien E, Parati G, Stergiou G, Asmar R, Beilin L, Bilo G, European Society of Hypertension Working Group on Blood Pressure Monitoring. European Society of Hypertension position paper on ambulatory blood pressure monitoring. *J Hypertens* 2013 Sep;31(9):1731-1768. [doi: 10.1097/HJH.0b013e328363e964] [Medline: 24029863]
36. [Declaración de Helsinki. Principios éticos para las investigaciones médicas en seres humanos]. Fortaleza, Brasil; 2013 Oct. URL: <http://www.wma.net/es/30publications/10policies/b3/index.html> [accessed 2017-12-04] [WebCite Cache ID 6vSfpALFh]
37. Topouchian J, Agnoletti D, Blacher J, Youssef A, Ibanez I, Khabouth J, et al. Validation of four automatic devices for self-measurement of blood pressure according to the international protocol of the European Society of Hypertension. *Vasc Health Risk Manag* 2011;7:709-717 [FREE Full text] [doi: 10.2147/VHRM.S27193] [Medline: 22174581]
38. Prieto L, Lamarca R, Casado A. [Assessment of the reliability of clinical findings: the intraclass correlation coefficient]. *Med Clin (Barc)* 1998 Feb 07;110(4):142-145. [Medline: 9541905]
39. Landis J, Koch G. The measurement of observer agreement for categorical data. *Biometrics* 1977 Mar;33(1):159-174. [Medline: 843571]
40. Portney L, Watkins M. *Foundations of Clinical Research: Applications to Practice*. 3rd edition. Upper Saddle River, NJ: Pearson-Prentice Hall; 2009.
41. Bland JM, Altman DG. Measurement error and correlation coefficients. *BMJ* 1996 Jul 06;313(7048):41-42 [FREE Full text] [Medline: 8664775]
42. Bland JM, Altman DG. Statistical methods for assessing agreement between two methods of clinical measurement. *Lancet* 1986 Feb 08;1(8476):307-310. [Medline: 2868172]
43. Topouchian J, Zelveyan P, Hakobyan Z, Melkonyan A, Asmar R. Validation of the QARDIO QARDIOARM upper arm blood pressure monitor, in oscillometry mode, for self measurement in persons fulfilling the population as described in this paper, according to the European Society of Hypertension International Protocol revision 2010. Dublin: dableEducational Trust; 2016. URL: <http://www.dableducational.org/Publications/2016/ESH-IP%202010%20Validation%20of%20QARDIO%20QARDIOARM.pdf> [WebCite Cache ID 6vMR7G66G]
44. Treskes RW, Wolterbeek R, van der Velde ET, Eindhoven DC, Schaliq MJ. Comparison of the diagnostic accuracy of four smartphone-compatible blood pressure monitors in post-myocardial infarction patients. *J Telemed Telecare* 2017 Jan 01:e1. [doi: 10.1177/1357633X17704092] [Medline: 28457182]
45. Bramlage P, Deutsch C, Krüger R, Wolf A, Müller P, Zwingers T, et al. Validation of the custo screen 400 ambulatory blood pressure-monitoring device according to the European Society of Hypertension International Protocol revision 2010. *Vasc Health Risk Manag* 2014;10:303-309 [FREE Full text] [doi: 10.2147/VHRM.S63602] [Medline: 24868162]
46. Benetti E, Fania C, Márquez HV, Palatini P. Validation of the Thermor BIOS BD215 device for home blood pressure measurement according to the European Society of Hypertension International Protocol revision 2010. *Blood Press Monit* 2014 Jun;19(3):176-179. [doi: 10.1097/MBP.0000000000000030] [Medline: 24589529]
47. Takahashi H, Yoshika M, Yokoi T. Validation of two automatic devices: Omron HEM-7252G-HP and Omron HEM-7251G for self-measurement of blood pressure according to the European Society of Hypertension International Protocol revision 2010. *Blood Press Monit* 2015 Oct;20(5):286-290. [doi: 10.1097/MBP.0000000000000127] [Medline: 25932887]
48. Wu N, Zhang X, Wang W, Zhang H. Validation of the Andon KD595 for clinical use and self-measurement according to the European Society of Hypertension International Protocol. *Blood Press Monit* 2016 Apr;21(2):124-127. [doi: 10.1097/MBP.0000000000000164] [Medline: 26544523]
49. Bing S, Chen K, Hou H, Zhang W, Li L, Wei J, et al. Validation of the Microlife BP A200 Comfort and W2 Slim automated blood pressure monitors in a general adult population according to the European Society of Hypertension and the ANSI/AAMI/ISO 81060-2: 2013 protocols. *Blood Press Monit* 2016 Apr;21(2):118-123. [doi: 10.1097/MBP.0000000000000169] [Medline: 26683381]

Abbreviations

BMI: body mass index

BP: blood pressure
DBP: diastolic blood pressure
ESH: European Society of Hypertension
HR: heart rate
ICC: intraclass correlation coefficient
SBP: systolic blood pressure
SEM: standard error of measurement

Edited by G Eysenbach; submitted 13.07.17; peer-reviewed by C Calvo Lobo, D López López, P Wark; comments to author 31.08.17; revised version received 03.10.17; accepted 30.10.17; published 15.12.17

Please cite as:

Mazoterias Pardo V, Losa Iglesias ME, López Chicharro J, Becerro de Bengoa Vallejo R
The QardioArm App in the Assessment of Blood Pressure and Heart Rate: Reliability and Validity Study
JMIR Mhealth Uhealth 2017;5(12):e198
URL: <http://mhealth.jmir.org/2017/12/e198/>
doi: [10.2196/mhealth.8458](https://doi.org/10.2196/mhealth.8458)
PMID: [29246880](https://pubmed.ncbi.nlm.nih.gov/29246880/)

©Victoria Mazoterias Pardo, Marta E Losa Iglesias, José López Chicharro, Ricardo Becerro de Bengoa Vallejo. Originally published in JMIR Mhealth and Uhealth (<http://mhealth.jmir.org>), 15.12.2017. This is an open-access article distributed under the terms of the Creative Commons Attribution License (<https://creativecommons.org/licenses/by/4.0/>), which permits unrestricted use, distribution, and reproduction in any medium, provided the original work, first published in JMIR mhealth and uhealth, is properly cited. The complete bibliographic information, a link to the original publication on <http://mhealth.jmir.org/>, as well as this copyright and license information must be included.